

兔抗 CUX1 多克隆抗体

中文名称：兔抗 CUX1 多克隆抗体

英文名称：Anti-CUX1 rabbit polyclonal antibody

别名：cut like homeobox 1 CDP; CUX; p75; CASP; CDP1; COY1; Clox; p100; p110; p200; CUTL1; GOLIM6; CDP/Cut; Cux/CDP; Nbla10317

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CUX1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the homeodomain family of DNA binding proteins. It may regulate gene expression, morphogenesis, and differentiation and it may also play a role in the cell cycle progression. Several alternatively spliced transcript variants encoding different isoforms have been identified.
Applications:	ELISA, IHC
Name of antibody:	CUX1
Immunogen:	Synthetic peptide of human CUX1
Full name:	cut like homeobox 1
Synonyms:	CDP; CUX; p75; CASP; CDP1; COY1; Clox; p100; p110; p

	200; CUTL1; GOLIM6; CDP/Cut; Cux/CDP; Nbla10317
SwissProt:	P39880
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human colorectal cancer
IHC Recommend dilution:	25-100



