

兔抗 DCP2 多克隆抗体

中文名称：兔抗 DCP2 多克隆抗体

英文名称：Anti-DCP2 rabbit polyclonal antibody

别名：decapping mRNA 2 NUDT20

抗原：DCP2

储存：冷冻（-20℃）

宿主：Rabbit

克隆类型：rabbit polyclonal

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

技术规格

Background:	The protein encoded by this gene is a key component of an mRNA-decapping complex required for degradation of mRNAs, both in normal mRNA turnover, and in nonsense-mediated mRNA decay (NMD). It removes the 7-methyl guanine cap structure from mRNA, prior to its degradation from the 5' end. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.
Applications:	ELISA, IHC
Name of antibody:	DCP2
Immunogen:	Fusion protein of human DCP2
Full name:	decapping mRNA 2

Synonyms:	NUDT20
SwissProt:	Q8IU60
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human esophagus cancer
IHC Recommend dilution:	25-100



