Recombinant Human DKK-1

Catalog Number: mly00701

Characteristics:

Source Chinese Hamster Ovary cell line

Description Human Dkk-1 (Thr32-His266)

Accession # O94907

Predicted molecular

mass

25.7kDa

Specification:

Formulation

Appearance White powder, Colorless clear liquid after reconstitution

Purity ≥95%, by SDS-PAGE (under reducing (R) & Non-reducing conditions,

visualized by Coomassie staining)

Endotoxin ≤10 EU/mg by the LAL method

Activity Measured by its ability to inhibit 50 ng/ml wnta-3a induced Topflash

reporter activity in CHO cells. The ED₅₀ is \leq 50ng/ml.

Lyophilized from a 0.22 µm-filtered solution containing PBS, 5% mannitol

and 0.01% Tween 80, pH 7.4

Size 50 μg/vial, 1 mg/vial

Handling and Storage:

Reconstitution It is recommended to redissolve in sterile deionized water.


Shipping The product is shipped with blue ice.

Storage & Stability 36 months at -20°C to -80°C in lyophilized state

6 months at -20°C to -80°C under sterile conditions after reconstitution 7-10 days at 2°C to 8°C under sterile conditions after reconstitution Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

Data:

Bioactivity



Measured by its ability to inhibit 50ng/ml Wnt-3a induced Top flash reporter activity in CHO cells.