

Recombinant Anti-Human CD52 mAb

Catalog Number : mIC02411

Characteristics:

Source	Chinese Hamster Ovary cell line
Description	Isotype: Human IgG1/Kappa
Predicted molecular mass	150 kDa

Specification:


Appearance	Colorless clear liquid
Purity	≥ 95%, by SDS-PAGE (under reducing (R) & Non-reducing conditions, visualized by Coomassie staining)
Endotoxin	≤1 EU/mg by the LAL method
Activity	Measure its binding ability with NK92 cells using flow cytometry.
Formulation	0.22 µm-filtered solution containing PBS, 0.1% Tween, pH 7.4
Size	1 mg/vial

Handling and Storage:

Shipping	The product is shipped with blue ice
Storage & Stability	36 months at -20°C to -80°C 7-10 days at 2°C to 8°C Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

Data:

Bioactivity



Recombinant Anti-Human CD52 mAb (Cat. No. mIC02411) binding with NK92 cells.