

www.mlbio.cn

mlbio

Anti-Dexamethasone Monoclonal Antibody

抗地塞米松单克隆抗体

Product name : Anti-Dexamethasone Monoclonal Antibody

Code : ml108808

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC₅₀ : 0.1ppb

CollaurationSensitivity : 0.3ppb

Protein Concentration : 7mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

Related products :

anti-Erythropoetin monoclonal antibody

anti-Calcitonin monoclonal antibody

anti-Cortisol monoclonal antibody

anti-Testosterone monoclonal antibody

anti-Prostaglandin monoclonal antibody

anti-Estrone monoclonal antibody

anti-Progesterone monoclonal antibody

anti-Estrone monoclonal antibody

anti-Estradiol monoclonal antibody

anti-Growth hormone monoclonal antibody

anti-Human chorionic gonadotropin monoclonal antibody

anti-uteinizing hormone/Prolactin monoclonal antibody

anti-Bisphenol A monoclonal antibody, anti-Thyroxin(e) monoclonal antibody

www.mlbio.cn



微信公众号



手机官网

企业信息:

上海酶联生物科技有限公司

办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼

电 话: 4008-898-798

传 真: 021-66980655

E-MAIL: 2881505699@qq.com