

www.mlbio.cn

mlbio

Anti-2-Nitrofurantoin metabolite Monoclonal Antibody

抗 2-硝基呋喃妥因代谢物单克隆抗体

Product name : Anti-2-Nitrofurantoin metabolite Monoclonal Antibody

Code : ml108859

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC₅₀ : 0.6ppb

ELISA Sensitivity : 0.1ppb

CollaaurumSensitivity : 5ppb

Protein Concentration : 6mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

Relatedproducts :

anti-2-Nitrofurazone metabolite monoclonal antibody

www.mlbio.cn

anti-2-Furaltadonemetabolite monoclonal antibody

anti-2-Furazolidone metabolite monoclonal antibody



微信公众号



手机官网

企业信息:

上海酶联生物科技有限公司

办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼

电 话: 4008-898-798

传 真: 021-66980655

E-MAIL: 2881505699@qq.com