



Anti- Cyclosporin A Monoclonal Antibody

抗环孢素 A 单克隆抗体

Product name : Anti- Cyclosporin A Monoclonal Antibody

Code : ml108911

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC₅₀ : 2ppb

ELISA Sensitivity : 0.2ppb

CollaaurumSensitivity : 10ppb

Protein Concentration : 7.3mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

Related products :

anti-Ampicillin monoclonal antibody

anti-Ceftiofur monoclonal antibody

anti-Cefquinome monoclonal antibody

anti-Cephalosporin monoclonal antibody

anti-Cephalexin monoclonal antibody

anti-Hygromycin B monoclonal antibody

anti-Bacitracin monoclonal antibody

anti-Polymyxin B monoclonal antibody

anti-Ampicillin Trihydrate monoclonal antibody

anti-metoclopramide monoclonal antibody

anti-Trimethoprim monoclonal antibody

anti-Polymyxin E monoclonal antibody

www.mlbio.cn



微信公众号



手机官网

企业信息:

上海酶联生物科技有限公司

办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼

电 话: 4008-898-798

传 真: 021-66980655

E-MAIL: 2881505699@qq.com