

www.mlbio.cn

mlbio

Diethylstilbestrol(Hexestrol Antigen)

己烯雌酚抗原

Product name : Diethylstilbestrol/Hexestrol antigen

Code : ml108990

Protein carrier : BALB/C mouse

Product type : BSA/OVA

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 30000

Collaurum titer : > 1 : 40, coupling ratio > 1 : 60

Protein concentration : 7.0mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 2 years, avoid repeated freezing and thawing

www.mlbio.cn



微信公众号



手机官网

企业信息:

上海酶联生物科技有限公司

办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼

电 话: 4008-898-798

传 真: 021-66980655

E-MAIL: 2881505699@qq.com