

www.mlbio.cn

mlbio

Polymyxin B Antigen

多粘菌素 B 抗原

Product name : Polymyxin B Antigen

Code : ml109584

Protein carrier : BALB/C mouse

Product type : BSA/OVA

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 30000

Collauration titer : > 1 : 40, coupling ratio > 1 : 60

Protein concentration : 8.0mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20°C

Validity : 2years, avoid repeated freezing and thawing

www.mlbio.cn



微信公众号



手机官网

企业信息:

上海酶联生物科技有限公司

办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼

电 话: 4008-898-798

传 真: 021-66980655

E-MAIL: 2881505699@qq.com