

www.mlbio.cn

mlbio

Anti-Hydrocortisone Monoclonal Antibody

抗氢化可的松单克隆抗体

Product name : Anti-Hydrocortisone Monoclonal Antibody

Code : ml300298

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC₅₀ : 3ppb

ELISA Sensitivity : 0.3ppb

CollaaurumSensitivity : 25ppb

Protein Concentration : 8mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

Related products :

anti-Tolfenamicacid monoclonal antibody

anti-Carprofen monoclonal antibody

anti-Furaltadone monoclonal antibody

anti-Metamizole Sodium Tablets monoclonal antibody

anti-Clonazepam monoclonal antibody

anti-Diazepam monoclonal antibody

anti-Antipyrinemonoclonal antibody

anti-Indometacin Tablets monoclonal antibody

www.mlbio.cn



微信公众号



手机官网

企业信息:

上海酶联生物科技有限公司

办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼

电 话: 4008-898-798

传 真: 021-66980655

E-MAIL: 2881505699@qq.com