



Anti-Ractopamine Monoclonal Antibody(RCT)

抗雷克托巴明单克隆抗体 (RCT)

Product name : Anti-Ractopamine Monoclonal Antibody

Code : ml300320

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC₅₀ : 0.5ppb

ELISA Sensitivity : 0.1ppb

CollaurumSensitivity : 3ppb

Protein Concentration : 6mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

Related products :

anti-Clenbuterol monoclonal antibody

www.mlbio.cn

anti-Salbutamol monoclonal antibody

anti-Chlorpromazine monoclonal antibody

anti-Cimaterol monoclonal antibody

anti-Cyproheptadine monoclonal antibody



微信公众号



手机官网

企业信息:

上海酶联生物科技有限公司

办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼

电 话: 4008-898-798

传 真: 021-66980655

E-MAIL: 2881505699@qq.com