



Anti-Aflatoxin M1 Monoclonal Antibody

抗黄曲霉毒素 M1 单克隆抗体

Product name : Anti-Aflatoxin M1 Monoclonal Antibody

Code : ml300336

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC₅₀ : 1ppb

CollaurationSensitivity : 0.5ppb

Protein Concentration : 8mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

Related products :

anti-T-2 toxin monoclonal antibody

anti-Aflatoxin B1 monoclonal antibody

anti-Total aflatoxin monoclonal antibody

anti-Ochratoxin A monoclonal antibody

anti-Zearalenone monoclonal antibody

anti-Fumonisin monoclonal antibody

anti-Deoxynivalenol monoclonal antibody

anti-Citrinin monoclonal antibody

anti-Sterigmatocystin monoclonal antibody

anti-Alternariol monoclonal antibody

anti-Ricin monoclonal antibody

anti-Microcystin LR monoclonal antibody

anti-Diphtheria toxin monoclonal antibody

anti-Cholera toxin monoclonal antibody

anti-Clostridium botulinum toxin A/B/D/E monoclonal antibody

anti-E. coli heat-labile enterotoxin A, B monoclonal antibody

www.mlbio.cn



微信公众号



手机官网

企业信息:

上海酶联生物科技有限公司

办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼

电 话: 4008-898-798

传 真: 021-66980655

E-MAIL: 2881505699@qq.com