

www.mlbio.cn

mlbio

Anti-6-Benzylaminopurine Monoclonal Antibody

抗 6 苯甲基氨基嘌呤单克隆抗体

Product name : Anti-6-Benzylaminopurine Monoclonal Antibody

Code : ml300363

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC₅₀ : 2-3ppb

ELISA Sensitivity: 0.3ppb

CollaurumSensitivity : 200ppb

Protein Concentration : 4.5mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

Related products :

anti-Imidacloprid monoclonal antibody

- anti-Carbendazim monoclonal antibody
- anti-Chlorothalonil monoclonal antibody
- anti-Methyl Parathion monoclonal antibody
- anti-Chlorpyrifos monoclonal antibody
- anti-Carbaryl monoclonal antibody
- anti-Aceto chlor monoclonal antibody
- anti-Class of Aceto chlor monoclonal antibody
- anti-Acetamiprid monoclonal antibody
- anti-Nitenpyram monoclonal antibody
- anti-Carbofuran monoclonal antibody
- anti-Aldicarb monoclonal antibody
- anti-Isoprocarb monoclonal antibody
- anti-Isoprocarb monoclonal antibody
- anti-4-Chlorophenoxyacetic acid monoclonal antibody
- anti-6-Benzylaminopurine monoclonal antibody
- anti-Paclobutrazol monoclonal antibody
- anti-Paraquat monoclonal antibody
- anti-Arermectin monoclonal antibody

www.mlbio.cn



微信公众号



手机官网

企业信息:

上海酶联生物科技有限公司

办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼

电话: 4008-898-798

传真: 021-66980655

E-MAIL: 2881505699@qq.com