

[www.mlbio.cn](http://www.mlbio.cn)

**mlbio**

## **Anti-Chlorothalonil Monoclonal Antibody**

### **抗氯噻酮单克隆抗体**

Product name : Anti-Chlorothalonil Monoclonal Antibody

Code : ml300425

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC<sub>50</sub> : 2-3ppb

ELISA Sensitivity: 0.3ppb

CollaurationSensitivity : 200ppb

Protein Concentration : 8mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

#### **Related products :**

anti-Imidacloprid monoclonal antibody

anti-Carbendazim monoclonal antibody

anti-Chlorothalonil monoclonal antibody

anti-Methyl Parathion monoclonal antibody

anti-Chlorpyrifos monoclonal antibody

anti-Carbaryl monoclonal antibody

anti-Aceto chlor monoclonal antibody

anti-Class of Aceto chlor monoclonal antibody

anti-Acetamiprid monoclonal antibody

anti-Nitenpyram monoclonal antibody

anti-Carbofuran monoclonal antibody

anti-Aldicarb monoclonal antibody

anti-Isoprocarb monoclonal antibody

anti-Isoprocarb monoclonal antibody

anti-4-Chlorophenoxyacetic acid monoclonal antibody

anti-6-Benzylaminopurine monoclonal antibody

anti-Paclobutrazol monoclonal antibody

anti-Paraquat monoclonal antibody

anti-Arermectin monoclonal antibody

[www.mlbio.cn](http://www.mlbio.cn)



微信公众号



手机官网

**企业信息:**

**上海酶联生物科技有限公司**

**办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼**

**电 话: 4008-898-798**

**传 真: 021-66980655**

**E-MAIL: 2881505699@qq.com**