

[www.mlbio.cn](http://www.mlbio.cn)

**mlbio**

## **Anti-Methyl Parathion Monoclonal Antibody**

### **抗甲基对硫磷单克隆抗体**

Product name : Anti-Methyl Parathion Monoclonal Antibody

Code : ml300432

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC<sub>50</sub> : 2-3ppb

ELISA Sensitivity: 0.3ppb

CollaurumSensitivity : 200ppb

Protein Concentration : 8mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

#### **Related products :**

anti-Imidacloprid monoclonal antibody

- anti-Carbendazim monoclonal antibody
- anti-Chlorothalonil monoclonal antibody
- anti-Methyl Parathion monoclonal antibody
- anti-Chlorpyrifos monoclonal antibody
- anti-Carbaryl monoclonal antibody
- anti-Aceto chlor monoclonal antibody
- anti-Class of Aceto chlor monoclonal antibody
- anti-Acetamiprid monoclonal antibody
- anti-Nitenpyram monoclonal antibody
- anti-Carbofuran monoclonal antibody
- anti-Aldicarb monoclonal antibody
- anti-Isoprocarb monoclonal antibody
- anti-Isoprocarb monoclonal antibody
- anti-4-Chlorophenoxyacetic acid monoclonal antibody
- anti-6-Benzylaminopurine monoclonal antibody
- anti-Paclobutrazol monoclonal antibody
- anti-Paraquat monoclonal antibody
- anti-Arermectin monoclonal antibody

[www.mlbio.cn](http://www.mlbio.cn)



微信公众号



手机官网

**企业信息:**

**上海酶联生物科技有限公司**

**办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼**

**电话: 4008-898-798**

**传真: 021-66980655**

**E-MAIL: 2881505699@qq.com**