



Anti-Paclobutrazol Monoclonal Antibody

抗帕洛司琼单克隆抗体

Product name : Anti-Paclobutrazol Monoclonal Antibody

Code : ml300434

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC₅₀ : 4-5ppb

ELISA Sensitivity: 0.5ppb

CollaurationSensitivity : 200ppb

Protein Concentration : 6mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

Related products :

anti-Imidacloprid monoclonal antibody

anti-Carbendazim monoclonal antibody

anti-Chlorothalonil monoclonal antibody

anti-Methyl Parathion monoclonal antibody

anti-Chlorpyrifos monoclonal antibody

anti-Carbaryl monoclonal antibody

anti-Acetoachlor monoclonal antibody

anti-Class of Acetoachlor monoclonal antibody

anti-Acetamiprid monoclonal antibody

anti-Nitenpyram monoclonal antibody

anti-Carbofuran monoclonal antibody

anti-Aldicarb monoclonal antibody

anti-Isoprocarb monoclonal antibody

anti-Isoprocarb monoclonal antibody

anti-4-Chlorophenoxyacetic acid monoclonal antibody

anti-6-Benzylaminopurine monoclonal antibody

anti-Paclobutrazol monoclonal antibody

anti-Paraquat monoclonal antibody

anti-Avermectin monoclonal antibody

www.mlbio.cn



微信公众号



手机官网

企业信息:

上海酶联生物科技有限公司

办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼

电 话: 4008-898-798

传 真: 021-66980655

E-MAIL: 2881505699@qq.com