

[www.mlbio.cn](http://www.mlbio.cn)

**mlbio**

## **Anti- Class of Albendazole Monoclonal Antibody**

### **抗阿苯达唑 I 类单克隆抗体**

Product name : Anti- Class of Albendazole Monoclonal Antibody

Code : ml300447

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC<sub>50</sub> : 10ppb

ELISA Sensitivity : 1ppb

CollaurationSensitivity : 500ppb

Protein Concentration : 8mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

#### **Related products :**

anti-Metronidazole monoclonal antibody

anti-4,4'- dinitrocarbanilide monoclonal antibody

anti- Natamycin monoclonal antibody

anti- Salinomycin monoclonal antibody

anti- Monensi monoclonal antibody

anti-Vancomycin monoclonal antibody

anti- Albendazole metabolite monoclonal antibody

anti-Maduramicin monoclonal antibody

anti- Diclazuril monoclonal antibody

[www.mlbio.cn](http://www.mlbio.cn)



微信公众号



手机官网

**企业信息:**

**上海酶联生物科技有限公司**

**办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼**

**电 话: 4008-898-798**

**传 真: 021-66980655**

**E-MAIL: 2881505699@qq.com**