

www.mlbio.cn

mlbio

Anti-Tylosin Monoclonal Antibody

抗托来霉素单克隆抗体

Product name : Anti-Tylosin Monoclonal Antibody

Code : ml300473

Source : BALB/C mouse

Product Type : IgG2a

Physical properties : Liquid or solid, avoid repeated freezing and thawing

ELISA titer : 0.005M PBS, > 1 : 400000

IC₅₀ : 0.5ppb

ELISA Sensitivity : 0.1ppb

CollaaurumSensitivity : 43ppb

Protein Concentration : 4.5mg/ml

Buffer : 0.01M PBS PH7.4

Storage : Below- 20 °C

Validity : 5 years, kept below 70 °C

Related products :

anti-Lincomycin monoclonal antibody

anti-Pirlimycin monoclonal antibody

anti-Erythromycin monoclonal antibody

anti-Avermectin monoclonal antibody

anti-Doramectin monoclonal antibody

anti-Ivermectin monoclonal antibody

anti-Tilmicosin monoclonal antibody

anti-Tiamulin monoclonal antibody

www.mlbio.cn



微信公众号



手机官网

企业信息:

上海酶联生物科技有限公司

办公地址: 上海市松江区云凯路66号临港科技绿洲T2栋16楼

电 话: 4008-898-798

传 真: 021-66980655

E-MAIL: 2881505699@qq.com