


CMTM3 Rabbit pAb

Host Species	Rabbit
Reactivity	Human, Monkey
Applications	WB, ELISA
MW	22kD (Observed)
Isotype	IgG
Storage	-15° C to -25° C/1 year (Do not lower than -25° C)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
Recommended Dilution Ratios	
WB 1:500-1:2000	
ELISA 1:40000	
Not yet tested in other applications.	
Basic Information	
Clonality	Polyclonal
Immunogen	The antiserum was produced against synthesized peptide derived from human CKLF3. AA range:1-50
Specificity	CMTM3 Polyclonal Antibody detects endogenous levels of CMTM3 protein.
Gene name	CMTM3
Protein Name	CKLF-like MARVEL transmembrane domain-containing protein 3


Organism	Gene ID	UniProt ID
Human	123920;	Q96MX0;
Mouse		Q99LJ5;

Cellular	Membrane;Multi-pass membrane protein.
Tissue specificity	Expressed in the leukocytes , placenta and testis.
Function	similarity:Belongs to the chemokine-like factor family. ,similarity:Contains 1MARVEL

domain. ,tissue speciﬁcity:Expressed in the leukocytes ,
placenta and testis. ,



Western Blot analysis of various cells using CMTM3 Polyclonal Antibody diluted at 1:1000



Western blot analysis of lysates from Jurkat and COS cells, using CKLF3 Antibody. The lane on the right is blocked with the synthesized peptide.